


## FACSVerse™ – Changing parameters in the Worklist

### During Experiment Acquisition in the Worklist

- Start acquisition.
- If you need to adjust PMT voltages, thresholds or modify gates, click **Stop Timer** to manually stop the acquisition. During preview mode, you can adjust threshold and PMT voltage and the system will automatically adjust spillover values for the current worklist only. You can also add or remove parameters before you run a worklist.

If you make changes, the Preview dialog opens. Select the third option if you want to apply these changes to all tubes, and then click OK.




- Click Resume to measure your sample.

*Attention:* In preview mode sample is used. Take care that there is enough sample left for measuring!

- ☑ If compensation is not correct and you need to adjust spillover values, click **Stop Timer** to manually stop the acquisition. During preview mode, you can adjust spillover values for the current worklist only.

If you make changes, the Preview dialog opens. Select the third option if you want to apply these changes to all tubes, and then click OK.



- ☑ Click Resume to measure your sample.

*Attention:* In preview mode sample is used. Take care that there is enough sample left for measuring.

*Remind:* These changes will not be saved in your assay they are only applied to the tubes in your worklist.