


|  |   |            |    |
|--|---|------------|----|
| <br>Proteomics-Facility<br>University of Konstanz | <b>Colloidal Coomassie<br/>         staining for<br/>         proteomic-samples</b> | SOP Nr.    | 00 |
|  |   | Pages      | 1  |
| Date   |   | 03-04-2013 |    |
| Version  |   | 1          |    |
| <b>SOP</b>   |   |            |    |

**Protocol is a Modification of:** Neuhoff, V., et al., *Improved staining of proteins in polyacrylamide gels including isoelectric focusing gels with clear background at nanogram sensitivity using Coomassie Brilliant Blue G-250 and R-250*. Electrophoresis, 1988. **9**(6): p. 255-62

## 1. Stock solutions:

### Solution A:

10% (w/v) ammonium sulphate, 2% (w/v) phosphoric acid in MilliQ water

### Solution B:

5% (w/v) Coomassie Brilliant Blue G-250 in MilliQ water

### Washing solution:

25% methanol in MilliQ water

## 2. Protocol

The gel is incubated for 5 min in water and then stained overnight in staining solution. The detailed protocol is described below:

- The staining solution is prepared by mixing 100 ml of the stock solution A with 2.5 ml stock solution B. Shake strongly this solution for 20 minutes and then add 25 ml pure methanol. Continue shaking the next 20-30 minutes.
- Place the gel in the freshly prepared colloidal Coomassie stain. Stain the gel overnight with gentle shaking.
- Wash the gel with 25 % methanol in MilliQ water for around 1 hour.
- Wash the gel with MilliQ water
- Repeat the last two steps until the protein bands are at the desired contrast against the background of the gel.